

EECS4312 Isolette Assignment

Michael Founk (michaelfounk@outlook.com)

Sean O'Brien (sean1025@my.yorku.ca)

December 10, 2018

Prism account used for submission: sean1025

©This document is not for public distribution. This document may only be used by EECS4312 students registered at York University. By downloading this document from the department, registered York students agree to keep this document (and all documents associated with assignments, projects or laboratories) private for their personal use, and may not communicate it to anyone else.

Students must obey York regulations on academic honesty requiring that students do the work of the Lab on their own, and not cheat by sharing with others or using and/or submitting the work of others. If you use *github* or similar repository for your work, the repository must be private. Placing your work in the public domain infringes on academic integrity. Github offers unlimited private repositories to students: <https://education.github.com/pack>.

Requirements Document:

Temperature control for an Isolette

Contents

1. System Overview	4
2. Goals	5
3. Context Diagram	6
4. Monitored Variables	7
5. Controlled Variables	8
6. Mode Diagram	9
7. R-Descriptions	10
8. E-descriptions	13
9. Abstract variables needed for the Function Table	15
10. Function Tables	16
10.1. Function Table for heat control: <i>c_hc</i>	16
10.2. Function Table for mode: <i>c_md</i>	17
10.3. Function Table for alarm: <i>c_al</i>	18
10.4. Function Table for temperature display: <i>c_td</i>	19
11. Use Case	20
A. Appendix	21

List of Figures

1. Isolette	4
2. Incubator Safety Problems [2, p98]	5
3. Context Diagram	6

List of Tables

1. Monitored Variables	7
2. Controlled Variables	8

1. System Overview

The System Under Development (SUD) is a computer controller for the thermostat of an Isolette.¹ An Isolette is an incubator for for an infant that provides controlled temperature, humidity and oxygen (Fig. 1). Isolettes are used extensively in Neonatal Intensive Care Units for the care of premature infants.

This requirements document is specifically for the control of temperature. The purpose of the Isolette computer controller is to maintain the air temperature of an Isolette within a desired range. It senses the current temperature of the Isolette and turns the heat source on and off to warm the air as needed. If the temperature falls too far below or rises too far above the desired temperature range, it activates an alarm to alert the nurse. The system allows the nurse to set the desired temperature range and to set the alarm temperature range outside the desired temperature range of which the alarm should be activated. This requirements documents follows the specification in [1] (Appendix A) except where noted.



Figure 1: Isolette

Many babies have died due to faulty incubators. There is thus a standard that manufacturers must satisfy. Modern incubators are equipped with alarms for air temperature, skin temperature, oxygen concentration and humidity. The alarms are both visual such

¹The image in Fig 1 is from: www.nufer-medical.ch.

as red warning lamps, and audio such as beep signals. Once measured values exceed permitted limits as well as when faults occur in sensors. For one such incident leading to death see “Medical Devices: Use and Safety” shown in Fig. 2.

CASE 6:2 Baby dies through overheating in incubator

An underdeveloped baby was being treated in an incubator with skin temperature control. When the baby was being washed, the skin sensor was removed and left hanging outside the incubator after the washing. Thus the sensor started measuring the room temperature (approx. 25°C). The control circuits therefore increased the heat to maximum level, and the temperature in the incubator rose to more than 45°C. The baby died.

For increased safety, incubators must be constructed with an extra control circuit that prevents overheating in case the skin sensor is misplaced. The incubator in question was indeed equipped with such a safety circuit, but the circuit was defective.

Figure 2: Incubator Safety Problems [2, p98]

2. Goals

The high-level goals (G) of the system are:

- G1—The Infant should be kept at a safe and comfortable temperature.
- G2—The Nurse should be warned if the Infant becomes too hot or too cold.
- G3—The cost of manufacturing the computer controller for the thermostat should be as low as possible.

3. Context Diagram

See Fig. A-1 in [1]. The System Under Description (SUD) is a computer *controller* to regulate the temperature of the Isolette. Everything else including the Operator Interface (described in [1]) is in the ecosystem (i.e. in the environment of the controller). The monitored variables and controlled variables for the controller are in Table 1 and Table 2, respectively. For clarity, simplicity and safety, there are some differences between the specifications in this document and the descriptions in [1].²




Figure 3: Context Diagram

²Documented in the write-up to this assignment: `assign1-spec.pdf`.

4. Monitored Variables

The monitored variables are a subset of those described in [1].³ There is a single status variable m_st that is *invalid* whenever any one of the operator inputs or temperature sensor are in a failed state. Otherwise types and ranges are as in [1].

Name	Type	Range	Units	Physical Interpretation
m_tm	\mathbb{R}	68 .. 105	°F	actual temperature of Isolette air temperature from sensor
m_dl	\mathbb{Z}	97 .. 99	°F	desired lower temperature set by operator
m_dh	\mathbb{Z}	98 .. 100	°F	desired higher temperature set by operator
m_al	\mathbb{Z}	93 .. 98	°F	lower alarm temperature set by operator
m_ah	\mathbb{Z}	99 .. 100	°F	higher alarm temperature set by operator
m_st	Enumerated	{valid, invalid}	State	status of sensor and operator settings
m_sw	Enumerated	{on, off}	State	switch set by operator

Table 1: Monitored Variables

³With some change of nomenclature. Monitored variables have an “m” prefix.

5. Controlled Variables

The controlled variables are a subset of those described in [1].⁴ In addition, there is a mode display c_md and a message display c_ms .⁵


Name	Type	Range	Units	Physical Interpretation
c_hc	Enumerated	{on, off}	State	heat control: command to turn heat source on or off
c_td	\mathbb{Z}	$\{0\} \cup \{68 \dots 105\}$	$^{\circ}\text{F}$	displayed temperature of Isolette (zero when Isolette is off)
c_al	Enumerated	{off, on}	State	sound alarm to call nurse
c_md	Enumerated	{off, init, normal, fail}	State	mode of Isolette operation (failed if $m_st = invalid$)
c_ms	Enumerated	{"ok", "er1", "er2", "er3"}	State	messages to display to nurse

Table 2: Controlled Variables

⁴With some change of nomenclature. Controlled variables have a “c” prefix.

⁵The mode “off” is added to that of Fig. A-4 in [1], and the mode transitions have been changed.

6. Mode Diagram



The system is in the off state before it is activated. Once it is turned on it enters the init state. In order for the system to enter the normal state, all three of the conditions must be met at the same time. When the system is in the normal state, and the sensor status is marked as invalid, the system moves into the fail state. If the sensor status is returned to valid, then the system moves from fail to normal. For all states, any permutation of variable values not including the ones specified beside an arrow result in the state remaining the same. Also the system can enter the Off state from any of the other three states; Init, Normal and Fail.

7. R-Descriptions

REQ1	The <i>controller</i> shall operate in one of four modes: <i>off</i> , <i>init</i> , <i>normal</i> and <i>fail</i> .	See mode diagram on page 6
------	--	----------------------------

REQ2	In the <i>normal</i> mode, the temperature controller shall maintain current temperature inside the Isolette within a set temperature range (the <i>desired</i> range).	The <i>desired</i> temperature range is $m_dl .. m_dh$. If the current temperature m_tm is outside this range, the controller shall turn the heater on or off via the controlled variable m_hc to maintain the desired state.
------	---	--

Rationale: The *desired temperature range* will be set by the nurse to the desired range based on the infant's weight and health. The controller shall maintain the current temperature within this range under normal operation.

The following relevant hazard was identified through the safety assessment process:

- **H1:** Prolonged exposure of Infant to unsafe heat or cold;
- *Classification:* catastrophic;
- *Probability:* $< 10^{-9}$ per hour of operation.

To ensure that probability of hazard H1 is 10^{-9} per hour of operation, the following derived safety requirement shall apply to the Isolette controller:

REQ3	<p>In <i>normal</i> mode, the controller shall activate an alarm whenever</p> <ul style="list-style-type: none">• the current temperature falls outside the <i>alarm</i> temperature range (either through temperature fluctuation or a change in the alarm range by an operator), or• a failure is signalled in any of the input devices (temperature sensor and operator settings).	<p>The alarm temperature range is $m_{al}..m_{ah}$. Monitored variable m_{st} shows “invalid” when any of the input signals fail.</p>
------	--	---

REQ4	<p>Once the alarm is activated, it becomes deactivated in one of two ways:</p> <ul style="list-style-type: none">• The nurse turns off the Isolette;• The alarm has lasted for 10 seconds, and after 10 seconds or more the alarm conditions are removed.	<p>Refer to the function table for alarm</p>
------	--	--

REQ5	In <i>normal</i> mode, if a sensor or heater malfunctions, the system transitions to the <i>fail</i> mode.	See mode diagram on page 6
------	--	----------------------------

Rationale: If a part of the system malfunctions the system is not operating correctly and must move to the *fail* mode.

REQ6	While the system is in the <i>fail</i> mode, the alarm is turned on.	See the function table for the alarm variable
------	--	---

Rationale: The operator must be alerted that the system is not currently functioning.

REQ7	While the system is in the <i>init</i> mode, the heater brings the temperature of the incubator into the desired range.	This is required for the system to move from "init" to "normal"
------	---	---

Rationale: Before the infant is put into the Isolette the temperature must meet the conditions supplied by the operator.

8. E-descriptions

ENV8	The current temperature received from the sensor is a real number in the range 68.0 to 105.0°F.	Refer to section 1
------	---	--------------------

Rationale: This is the specified range of operation of the Isolette. The lower end of this range is useful for monitoring an Isolette that is warming to the Desired Temperature Range. The upper end is set to be greater than the maximum Upper Alarm Temperature.

ENV9	The desired and alarm temperatures received from the operator are all in increments of 1°F.	Refer to the table of monitored variables
------	---	---

Rationale: Marketing studies have shown that customers prefer to set temperatures in 1 degree increments. A resolution 1°F is sufficient to be consistent with the functional and performance requirements specified in the rest of the document.

ENV10	An infant will only be placed in the Isolette in the normal mode	Refer to section 1
-------	--	--------------------

Rationale: The normal mode is the only state that guarantees the conditions that keep the infant safe.

ENV11	The lower alarm temperature will always be less than or equal to the lower desired temperature.	Refer to the monitored/controlled variables tables.
-------	---	---

Rationale: If the lower alarm temperature is greater than the lower desired temperature then the alarm can go off while the current temperature is still in the desired temperature range.

ENV12	The lower desired temperature will always be less than or equal to the upper desired temperature.	Refer to the monitored/controlled variables tables.
-------	---	---

Rationale: If the lower desired temperature is greater than the upper desired temperature it is unclear if the heat source should be on or off.

9. Abstract variables needed for the Function Table

Input Conditions	valid_temp(i)
$m_dl(i) < m_dh(i) \text{ AND } m_dl(i) \leq m_tm(i) \text{ AND } m_tm(i) \leq m_dh(i)$	TRUE
$\text{NOT}(m_dl(i) < m_dh(i) \text{ AND } m_dl(i) \leq m_tm(i) \text{ AND } m_tm(i) \leq m_dh(i))$	FALSE

Input Conditions	valid_desired(i)
$m_al(i) < m_dl(i) \text{ AND } m_dl(i) < m_dh(i) \text{ AND } m_dh(i) < m_ah(i)$	TRUE
$\text{NOT}(m_al(i) < m_dl(i) \text{ AND } m_dl(i) < m_dh(i) \text{ AND } m_dh(i) < m_ah(i))$	FALSE

10. Function Tables

Starting on the next page, provide one function table for each control variable (in Table 2). Each control variable should have its own sub-section heading and its own page.

10.1. Function Table for heat control: c_{hc}

Input Conditions					c_hc(i)
i = 0					off
i > 0	c_md(i-1) = off				off
	c_cm(i-1) = init	m_st(i) = invalid			off
		m_st(i) = valid	valid_desired(i)	m_tm(i) > m_dh(i)	off
				m_tm(i) < m_dl(i)	on
			NOT(valid_desired(i))		
	c_md(i-1) = normal	valid_desired(i)	m_tm(i) > m_dh(i)		off
			m_tm(i) < m_dl(i)		on
		NOT(valid_desired(i))			off
c_md(i-1) = fail				off	

10.2. Function Table for mode: c_md

Input Conditions				c_md(i)
i = 0				off
i > 0	m_sw(i-1) = off			off
	m_sw(i-1) = on	c_md(i-1) = off		init
		c_md(i-1) = init	m_st(i) = valid AND valid_desired(i) and valid_temp(i)	normal
			NOT(m_st(i) = valid AND valid_desired(i) and valid_temp(i))	c_md(i-1)
		c_md(i-1) = normal	m_st(i) = valid	c_md(i-1)
			m_st(i) = invalid	fail
		c_md(i-1) = fail	m_st(i) = valid	normal
			m_st(i) = invalid	c_md(i-1)

10.3. Function Table for alarm: c_al

Input Conditions					c_al
i = 0					off
i > 0	c_md(i-1) = off				off
	c_md(i-1) = init				off
	c_md(i-1) = normal	valid_desired(i) AND valid_temp(i)	c_al(i-1) = on	held_for(c_al,10)(i-1)=TRUE	off
				held_for(c_al,10)(i-1)=FALSE	on
			c_al(i-1) = off		off
		NOT(valid_desired(i)) OR NOT (valid_temp(i))			on
c_md(i-1) = fail				on	

10.4. Function Table for temperature display: c_td

Input Conditions		c_td
$i = 0$		0
$i > 0$	$c_md(i-1) = \text{normal}$	$m_tm(i)$
	$\text{NOT}(c_md(i-1) = \text{normal})$	0

11. Use Case

Related System Goals: G1

Primary Actor: Nurse

Precondition: - The Isolette is in normal mode with an infant inside - The alarm is off

Postcondition:

- The Isolette is in normal mode with new desired and alarm temperature ranges
- The alarm is off
- The infant is inside the Isolette

Main Success Scenario:

- The infant's healthy range for temperature must be adjusted
- The nurse removes the infant from the Isolette
- The nurse turns off the Isolette
- The nurse turns on the Isolette, and specifies the desired ranges
- The Isolette sets the temperature to the desired range
- The nurse places the infant in the Isolette

Alternate Course 1:

- Nurse removes infant from Isolette
- Nurse turns off the Isolette
- Nurse initializes another Isolette
- Nurse places infant in new Isolette when temperature reaches desired range

References

- [1] US FAA. Requirements Engineering Management Handbook. Technical Report DOT/FAA/AR-08/32, U.S. Department of Transportation Federal Aviation Administration, June 2009.
- [2] Bertil Jacobson and Alan Murray. *Medical Devices: Use and Safety*. Elsevier, 2007.

A. Appendix

```

Alarm    [ delta: posreal ]
          : THEORY

BEGIN

    importing Time[delta]
    importing Mode[delta]
    i: VAR DTIME

    STATE: NONEMPTY_TYPE = {valid, invalid}
    LOW_DESIRE: TYPE+ = subrange(97, 99)
    UPPER_DESIRE: TYPE+ = subrange(98, 100)
    LOW_ALARM: TYPE+ = subrange(93, 98)
    UPPER_ALARM: TYPE+ = subrange(99, 103)
    SWITCH: NONEMPTY_TYPE = {on, off}
    TEMP_RANGE: NONEMPTY_TYPE = {r: real | (68.0 <= r AND
    r <= 105.0) OR r = 0}
    MODE: NONEMPTY_TYPE = {off, init, normal, fail}

    m_st : [DTIME -> STATE]
    m_dl : [DTIME -> LOW_DESIRE]
    m_dh : [DTIME -> UPPER_DESIRE]
    m_al : [DTIME -> LOW_ALARM]
    m_ah : [DTIME -> UPPER_ALARM]
    m_sw : [DTIME -> SWITCH]
    m_tm : [DTIME -> TEMP_RANGE]

    c_al: [DTIME -> SWITCH]

    convert_tf(i): bool =
        COND
        c_al(i) = on -> TRUE,
        c_al(i) = off -> FALSE
        ENDCOND

    c_al_ft(i): bool =
        COND
        i = 0 -> c_al(i) = off,
        i > 0 ->

```

```

        COND
        c_md(i-1) = off -> c_al(i) = off,
        c_md(i-1) = init -> c_al(i) = off,
        c_md(i-1) = normal ->
            COND
            valid_desired(i) ->
                COND
                valid_temp(i) -> c_al(i)
                = off,
                NOT (valid_temp(i)) ->
                c_al(i) = on
                ENDCOND,
            NOT (valid_desired(i)) ->
            c_al(i) = on
            ENDCOND,
        c_md(i-1) = fail ->
            COND
            c_al(i-1) = on ->
                COND
                held_for(convert_tf,10) (i-1)
                = TRUE -> c_al(i) = off,
                held_for(convert_tf,10) (i-1)
                = FALSE -> c_al(i) = on
                ENDCOND,
            c_al(i-1) = off -> c_al(i) = on
            ENDCOND
        ENDCOND
    ENDCOND

END Alarm

```

```

HeatControl [ delta: posreal ]
: THEORY

```

```

BEGIN

```

```

    importing Time[delta]
    importing Mode[delta]
    i: VAR DTIME

```

```

c_hc : [DTIME -> SWITCH]

c_hc_func(i): bool =
    COND
    i = 0 -> c_hc(i) = off,
    i > 0 ->
        COND
        c_md(i-1) = off -> c_hc(i) = off,
        c_md(i-1) = init ->
            COND
            m_st(i) = invalid -> c_hc(i) = off,
            m_st(i) = valid ->
                COND
                NOT(valid_desired(i)) ->
                    c_hc(i) = off
                valid_desired(i) ->
                    COND
                    m_tm(i) > m_dh(i)
                    -> c_hc(i) = off,
                    m_tm(i) < m_dl(i)
                    -> c_hc(i) = on
                ENDCOND
            ENDCOND
        c_md(i-1) = normal ->
            COND
            valid_desired(i) ->
                COND
                m_tm(i) > m_dh(i) ->
                    c_hc(i) = off,
                m_tm(i) < m_dl(i) ->
                    c_hc(i) = on
                ENDCOND
            NOT(valid_desired(i)) ->
                c_hc(i) = off
            ENDCOND
        c_md(i-1) = fail -> c_hc(i) = off
    ENDCOND
END HeatControl

```

```

TemperatureDisplay    [ delta: posreal ]

```

```

                                : THEORY

BEGIN

    importing Time[delta]
    importing Mode[delta]
    i: VAR DTIME

    c_td : [DTIME -> real]

    c_td_func(i): bool =
        COND
        i = 0 -> 0
        i > 0 ->
            COND
            c_md(i-1) = off -> c_td(i) = 0,
            c_md(i-1) = init -> c_td(i) = 0,
            c_md(i-1) = normal -> c_td(i) = m_tm(i),
            c_md(i-1) = fail -> c_td(i) = 0
            ENDCOND
        ENDCOND

END TemperatureMeter
```